

## Vector® TMB Peroxidase Substrate

**Cat. No.:** SK-4400

**Storage:** 2-8 °C

### DESCRIPTION

This kit contains all of the reagents necessary to prepare the working solution. Vector TMB Substrate (3, 3', 5, 5'-tetramethylbenzidine) produces a blue reaction product in the presence of peroxidase (HRP) enzyme. On tissue sections, TMB is more sensitive than other peroxidase substrates but the reaction product has a more diffuse appearance. By omitting the Stabilizing Solution from the TMB kit the working solution will produce a soluble blue product useful for ELISA.

### COMPONENTS

<u>Product Name</u>	<u>Volume</u>
Vector TMB Reagent 1	6 ml
Vector TMB Reagent 2	6 ml
Vector TMB Reagent 3	6 ml
Vector TMB Reagent 4	6 ml

### STORAGE:

- Store reagents in original bottles at 2-8 °C
- Avoid storing reagents or working solution in strong direct light


### PREPARATION OF SUBSTRATE WORKING SOLUTION:

- To 5.0 ml of distilled water
- Add 2 drops (≈ 72 μl<sup>†</sup>) of Vector TMB Reagent 1
- Add 3 drops (≈ 81 μl<sup>†</sup>) of Vector TMB Reagent 2
- Add 2 drops (≈ 80 μl<sup>†</sup>) of Vector TMB Reagent 3
- Add 2 drops (≈ 80 μl<sup>†</sup>) of the Vector TMB Reagent 4

-Mix well before use. Use immediately.

<sup>†</sup>Drop volumes differ due to solvent compositions.

**IMPORTANT:** Little is known about the toxicity and carcinogenicity of the substrate kit components. Appropriate care should be exercised when using this reagent including gloves, eye protection, lab coats, and good laboratory procedures. Dispose in accordance with local regulations.



*Tumor: Cytokeratin (s) detected with VECTASTAIN® Elite® ABC Kit and Vector TMB Substrate (blue). No counterstain.*

### INSTRUCTIONS FOR USE

#### For Tissues or Cells

After incubation with a peroxidase (HRP) detection system, rinse sections in buffer. Incubate with the substrate working solution for 5-15 minutes. Optimal development times should be determined by the investigator.

Wash slides with buffer for 2-3 minutes. Rinse in water and counterstain if desired (see chart on reverse). Excessive wash times after substrate reaction can cause decreased sensitivity.

**For permanent, non-aqueous mounting:** Dehydrate, clear and coverslip using a non-aqueous mounting media, such as VectaMount® Mounting Medium (H-5000). Do not mount with aqueous mounting media.

#### Use in ELISA

To use the Vector TMB Substrate for ELISA, do not add the Stabilizing Solution. Measure the absorbance of the soluble blue product at 650 nm. To increase absorbance of the blue Vector TMB product two to four-fold, stop the reaction by adding 50 μl of 1N sulfuric acid per microtiter plate well. Measure the absorbance of the yellow reaction product at 450 nm.

### NOTES

We recommend using glass-distilled water in the preparation of the substrate buffer. Deionized water may contain inhibitors of the peroxidase reaction.

Detailed product listings, specifications, protocols and additional information is available on our website: [vectorlabs.com](http://vectorlabs.com)